



Use of SNS Ventilators in the Pediatric Patient (LP-10, LTV-1200 and Uni-Vent Eagle)

Concepts of Rapid or Volume Based Pediatric Ventilation

- If volume ventilating, start at 10 ml/kg (unless protective lung strategy ventilation required)
- Volume lost to circuit must be replaced unless measurements taken at “wee”
- Set I-time generally between 0.7 – 1.0 sec.

NORMAL RESPIRATORY RATES

- Infant 30-60
- Toddler 24-40
- Preschooler 22-34
- School-age child 18-30
- Adolescent 12-16

ASSESSMENT

- Chest rise
- Breath sounds
- Respiratory Rate
- Work of Breathing
- Pressures required to deliver volume
- ABG/TCM/SaO₂



UNI-VENT® Eagle™ Ventilation System “Quick Set-Up”

- Mode** ▶ SIMV, Assist Control (A/C), CPAP
- Volume** ▶ 10 ml/kg (displayed on LCD above control)
- Breath rate** ▶ Set age appropriate (dial is sensitive, I-150 bpm)

- Inspiratory Time** ▶ (0.7-1.0 Combination of Inspiratory time and I:E ratio is displayed on LCD (I:E ratio default 1:1 preset)
- FiO₂** ▶ Dial desired FiO₂ (21% to 100%) value displayed on LCD
- PEEP** ▶ Pushbutton switch: each push = 1 cwp
- Alarms** ▶ Set based on average Peak Inspiratory Pressure
- Low Alarm Limit** ▶ 5 cwp below spontaneous Peak Inspiratory Pressure
- High Alarm Limit** ▶ 10 cwp above Mechanical Breaths Peak Inspiratory Pressure
- Battery Life** ▶ 3 hours maximum using internal compressor; 12 hours using external gas source



Use of SNS Ventilators in the Pediatric Patient (LP-10, LTV-1200 and Uni-Vent Eagle)



LP-10 Ventilator “Quick Set Up”

- | | |
|---|--|
| Mode | ▶ SIMV or Assist Control |
| Volume | ▶ 10ml/kg (measured by spirometer) |
| Breath rate | ▶ Set age appropriate |
| Inspiratory Time | ▶ (0.7-1.0 second) |
| FiO₂ | ▶ Capable of delivering 100% |
| Back of machine- O₂ enrichment kit | ▶ Set liter flow (not > 10 lpm)
▶ Analyze FiO ₂ |
| Front of machine- Bleed O₂ into circuit | ▶ Highest FiO ₂ approximately 40%
▶ Adjust liter flow and analyze |
| PEEP | ▶ External PEEP valve located on circuit exhalation valve; dial in desired PEEP pressure |
| Alarms | ▶ Set based on average Peak Inspiratory Pressure |
| Low Alarm Limit | ▶ 5 cwp below Peak Inspiratory Pressure |
| High Alarm Limit | ▶ 10 cwp above Mechanical Breaths Peak Inspiratory Pressure |
| Battery Life | ▶ Internal Battery 30 minutes to 1 hour; External battery 10 hours |



LTV-1200 “Quick Set Up”

- | | |
|------------------------|---|
| Preuse | ▶ Vent Op/Leak test/New Patient/Patient Size |
| Mode | ▶ SIMV, Assist Control (A/C), Pressure Support-CPAP, NPPV |
| Volume | ▶ 10 ml/kg |
| Breath Rate | ▶ Set age appropriate |
| PEEP | ▶ Set on Machine |
| FiO₂ | ▶ High pressure source – set oxygen on vent (capable of delivering 100%)
▶ Low pressure source – adjust liter flow and analyze |
| Alarms | ▶ Low Pressure – Set 5 cwp below average spontaneous Peak Inspiratory Pressure
▶ High Pressure – Set 10 cwp below average Peak Inspiratory Pressure
▶ Low Min. Vol. – Set 1 lpm below average minute volume |
| Ext. Features | ▶ Safe to use defaults provided appropriate size patient selected |
| Battery Life | ▶ Internal Battery 1 hour; Small external battery 3 hours; Large external battery 9 hours |